

---

# TruScanEEGpy Documentation

*Release 0.0.5*

**Dominique Makowski**

Oct 23, 2021



---

## Contents

---

<b>1 TruScanEEGpy</b>	<b>1</b>
1.1 Installation . . . . .	1
1.2 Documentation . . . . .	1
<b>2 Installation</b>	<b>3</b>
2.1 Stable release . . . . .	3
2.2 From sources . . . . .	3
<b>3 Functions</b>	<b>5</b>
<b>4 Indices and tables</b>	<b>7</b>
<b>Python Module Index</b>	<b>9</b>
<b>Index</b>	<b>11</b>



# CHAPTER 1

---

TruScanEEGpy

---

**Utility functions to work with Deymed's TruScan EEG system.**

## 1.1 Installation

To install TruScanEEGpy, run this command in your terminal:

```
$ pip install TruScanEEGpy
```

## 1.2 Documentation

```
> import TruScanEEGpy
> import mne
>
> layout = TruScanEEGpy.layout_128('10-5')
> montage = TruScanEEGpy.montage_mne_128(layout)
> mne.viz.plot_montage(montage)
```

docs/img/montage\_128.png

# CHAPTER 2

---

## Installation

---

### 2.1 Stable release

To install TruScanEEGpy, run this command in your terminal:

```
$ pip install TruScanEEGpy
```

This is the preferred method to install TruScanEEGpy, as it will always install the most recent stable release. It is also possible to install it directly from github:

```
$ pip install https://github.com/neuropsychology/TruScanEEGpy/zipball/master
```

If you don't have `pip` installed, this Python installation guide can guide you through the process.

### 2.2 From sources

The sources for TruScanEEGpy can be downloaded from the [Github repo](#).

You can either clone the public repository:

```
$ git clone git://github.com/neuropsychology/TruScanEEGpy
```

Or download the [tarball](#):

```
$ curl -OL https://github.com/neuropsychology/TruScanEEGpy/tarball/master
```

Once you have a copy of the source, you can install it with:

```
$ python setup.py install
```



# CHAPTER 3

---

## Functions

---

This part of the documentation details the complete TruScanEEGpy for python API.

Submodule.

**layout\_128 (names='index')**

Create layout file

**Parameters** **names** (*str*) – Can be ‘index’ or ‘10-5’ for the ‘mak-128’ (an adapted version of 10-5) layout.

### Examples

```
>>> import TruScanEEGpy
>>>
>>> layout = TruScanEEGpy.layout_128()
>>> layout = TruScanEEGpy.layout_128('10-5')
```

Submodule.

**montage\_mne\_128 (layout)**

Create MNE montage

### Examples

```
>>> import TruScanEEGpy
>>> import mne
>>>
>>> layout = TruScanEEGpy.layout_128('10-5')
>>> montage = TruScanEEGpy.montage_mne_128(layout)
>>> mne.viz.plot_montage(montage)
>>> mne.viz.plot_montage(montage, kind = "3d")
>>>
```

(continues on next page)

(continued from previous page)

```
>>> # Compare to other standard montages
>>> montage = mne.channels.make_standard_montage('standard_1005')
>>> montage.get_positions()
>>> montage.get_positions()['ch_pos']["TP9"]
>>> mne.viz.plot_montage(montage)
>>>
>>> montage = mne.channels.make_standard_montage('biosemi128')
>>> montage.get_positions()
>>> mne.viz.plot_montage(montage)
>>>
>>> montage = mne.channels.make_standard_montage('GSN-HydroCel-128')
>>> montage.get_positions()
>>> mne.viz.plot_montage(montage)
>>>
>>> montage = mne.channels.make_standard_montage('easycap-M1')
>>> mne.viz.plot_montage(montage)
```

# CHAPTER 4

---

## Indices and tables

---

- genindex
- modindex
- search



---

## Python Module Index

---

t

TruScanEEGpy.layout, 5  
TruScanEEGpy.montage, 5



---

## Index

---

### L

`layout_128()` (*in module* `TruScanEEGpy.layout`), 5

### M

`montage_mne_128()` (*in module* `TruScanEEGpy.montage`), 5

### T

`TruScanEEGpy.layout` (*module*), 5

`TruScanEEGpy.montage` (*module*), 5